

# Getting started with Coraye

After creating your account and activating your license, we will discuss the chapters on installing Coraye, installing plugins and presenting the software interface to you.

---

## Installation on Mac OS

64 bits compatible, minimum version Mac OS 10.10 (Yosemite).

For more information, see chapter: [Installation in Mac OS](#)

---

## Windows installation

x64 compatible (x32 version on request), minimum version Windows 7.

For more information, see the chapter: [Installation on Windows](#)

---

## Installing the Windows driver for your spectrophotometer

The drivers for Windows are supplied with the plugings. Regarding the Mac OS version, there is no driver to install.

For more information, see chapter: [Installing the Windows driver for your spectrophotometer](#)

---

## Installation and management of plugins

Don't forget to install the plug-in corresponding to your spectrophotometer. The management of plugins is done directly via the Coraye software.

---

To find out more, see the chapter: [Installing and managing plugins](#)

---

## Interface presentation

For more information, see chapter: [Interface presentation](#)

---

Revision #2

Created 6 September 2021 07:21:43 by Lionel WETTEREN

Updated 6 October 2021 14:49:40 by Lionel WETTEREN